

## Validation-, routine- and Process control sets

Due to the different requirements for the process control, the routine control and validation of processes, there are different sets at ebro.

As an example, applications for the different sets are clearly displayed in the table. So, it is possible to choose the set suitable for your application at a glance.

<b>SL 1011</b>	<b>SL 1111</b>	<b>SL 1201</b>	<b>SL 1521</b>	<b>SL 1621</b>	<b>SL 1221</b>	<b>SL 2001</b>	<b>SL 3001</b>	<b>SL 3101</b>	<b>SL 3301</b>	<b>SL 4011</b>	<b>SL 4102</b>	<b>SL 4121</b>
1250-1011	1250-1111	1250-1201	1250-1521	1250-1621	1250-1221	1250-2001	1250-3001	1250-3101	1250-3301	1250-4011	1250-4102	1250-4121
<b>Routine control in</b>												
Bedpan washer	✓											
Washer disinfecter		✓	✓		✓	✓						
Endoscope washer disinfecter		✓	✓			✓						
Sterilizer < 60 l		✓	✓	BD *	BD */ ✓	✓						
Sterilizer > 60 l		✓	✓	BD *	BD */ ✓	✓		✓		✓		
<b>Validation of Processes in</b>												
Washer disinfecter							✓			✓		
Endoscope washer disinfecter						✓	✓			✓		
Sterilizer < 60 l							✓	✓	✓	✓		
Sterilizer > 60 l									✓			
DAC Universal										✓		
<b>Process control in</b>												
Pasteurization cans										✓		
Pasteurization bottles										✓	✓	

\* Bowie&Dick-Test